

Semper
guard

force™ black



Serious Protection for Serious Work



Medical



Pesticides



First
Responder



Automotive



Beauty
& Tattoo



Manufacturing

*semperguard is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





force™ black

Semperguard® force™ black exam nitrile gloves are made from high-quality nitrile for durability and puncture resistance. Textured fingers provide a secure grip, ensuring users can handle tools and materials confidently, even in challenging conditions. These gloves are well-suited for industries like automotive, industrial work, tattoo artistry, and pesticides, where protection from chemicals, oils, and other contaminants is crucial, while still maintaining dexterity and comfort.

- Sleek and professional appearance.**
 Hides stains caused from grease, dirt, dyes, and ink
- ESD tested and approved**
 Safe for applications handling sensitive electronics according to ANSI/ESD STM 15.1
- Pesticide protection**
 Approved for handling diluted pesticides
- Tested against 36 chemotherapy and hazardous drugs, fentanyl, and gastric acid**
 Meets USP 800 requirements for PPE compliance for chemotherapy
- Safe for food contact applications**
 Meets FDA 21 CFR Part 177 requirements for food handling
- NFPA tested**
 2025 Certified: NFPA 1950 (1999)

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Black	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 5.1mils	Palm 3.5mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	572% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori-BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, ISO 18889, ASTM D6978, NFPA 1950	

Packaging

DIMENSIONS

Dispenser: 245L x 123W x 54H mm
Carton: 280L x 254W x 255H mm

Reorder numbers / Glove Counts*

	BKNF102	Small	100 gloves/bx	10 bxs/cs
	BKNF103	Medium	100 gloves/bx	10 bxs/cs
	BKNF104	Large	100 gloves/bx	10 bxs/cs
	BKNF105	X-Large	100 gloves/bx	10 bxs/cs
	BKNF106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



avg. 10° Q under worst-case test conditions.

SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks. For USP <800> chemotherapy use, change gloves every 20-30 minutes.

Licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at <https://gotgloves.com/products/>. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global