

Semper
guard

sapphire XP™ white



**Clean Protection,
Precise Control**



Beauty &
Tattoo



General
Purpose



Food &
Beverage



Non-Medical
Lab



Janitorial
& Sanitation



Safety &
Protection

*semperguard is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





sapphire XP™ white

Semperguard® sapphire XP™ white disposable nitrile gloves deliver dependable protection and exceptional grip, making them an excellent choice for tasks that demand both dexterity and safety. The unique, white color is suitable for industries such as janitorial-sanitation, safety, laboratory/pharma, and food services, where hygiene, cleanliness, and visibility are essential. Their superior fit and feel ensure maximum tactile sensitivity, allowing one to handle delicate tasks with precision and comfort.

- Clean white color**
 Supports a hygienic appearance and helps improve visibility of contamination.
- Comfortable fit**
 Soft nitrile formulation supports dexterity and tactile sensitivity.
- ESD tested and approved**
 Safe for applications handling sensitive electronics according to ANSI/ESD STM 15.1.
- Safe for food contact applications**
 Meets FDA 21 CFR Part 177 requirements for food handling.

Technical data

TYPE	Non-sterile disposable protective glove	
GRADE	General Purpose	
COLOR	White	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.9mils Palm 2.8mils	based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	575% (avg.)	
BARRIER PERFORMANCE	AQL < 4.0 (pinholes)	
SOCIAL COMPLIANCE	amfori-BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	EN ISO 21420	

Packaging

DIMENSIONS	Dispenser: 245L x 125W x 42H mm Carton: 260L x 223W x 255H mm
------------	--

Reorder numbers / Glove Counts *

	GSWNF102 Small	100 gloves/bx	10 bxs/cs
	GSWNF103 Medium	100 gloves/bx	10 bxs/cs
	GSWNF104 Large	100 gloves/bx	10 bxs/cs
	GSWNF105 X-Large	100 gloves/bx	10 bxs/cs

Certified Claims



avg. 10° D under worst-case test conditions.



SAMPLES: Available by request - ask your sales rep!
 *glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global