



# latex edition IC™



## Comfort and Strength



Ambulatory Surgical Centers



Dental



Veterinary



Long Term Care



Primary Care



Non-medical Lab

sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd.

CONTACT US! [gloveinfo@harpsglobal.com](mailto:gloveinfo@harpsglobal.com) | [www.gotgloves.com](http://www.gotgloves.com)

HARPS USA, Inc.  
3001 Executive Drive, Suite 330  
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.  
1 Harbourfront Avenue, #14-07 Keppel Bay Tower  
Singapore 098632, Singapore





# latex edition IC™

Sempermed® latex edition IC™ exam latex gloves are designed to provide the perfect combination of protection, tactile sensitivity and value. The advanced manufacturing technology and synthetic polymer inner coating creates a glove that is strong, flexible and goes on with ease.

- Specially formulated with **comfort and protection** in mind.
- Synthetic polymer inner coating allows for **easy donning**.
- The micro-textured surface finish offers a **strong grip** in wet or dry environments.

## Technical data






<b>TYPE</b>	Non-sterile latex examination and disposable protective glove	
<b>GRADE</b>	Medical	
<b>COLOR</b>	Natural	
<b>MATERIAL</b>	Made with natural rubber latex	
<b>GLOVE SURFACE</b>	Textured	
<b>COUNTRY OF ORIGIN</b>	Malaysia	
<b>SHelf LIFE</b>	5 years in original packaging if stored according to product specification	
<b>SIZE</b>	Length: 230mm (min) Width: 94mm +/- 5mm (based on size M)	
<b>WALL THICKNESS</b>	Finger 5.5mils	Palm 4.7mils
	based on average production value measured in single layer	
<b>TENSILE STRENGTH</b> <i>(before aging)</i>	20MPa (avg.)	
<b>ELONGATION</b> <i>(before aging)</i>	709% (avg.)	
<b>BARRIER PERFORMANCE</b>	AQL < 1.5 (pinholes)	
<b>SOCIAL COMPLIANCE</b>	amfori BSCI, EcoVadis	
<b>TESTING &amp; CONFORMITY TO STANDARDS</b>	ASTM D3578, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5	

## Packaging

### DIMENSIONS

Dispenser: 240L x 122W x 65H mm  
Carton: 340L x 252W x 250H mm

## Reorder numbers / Glove Counts\*

	<b>SCLT101</b>	X-Small	100 gloves/bx	10 bxs/cs
	<b>SCLT102</b>	Small	100 gloves/bx	10 bxs/cs
	<b>SCLT103</b>	Medium	100 gloves/bx	10 bxs/cs
	<b>SCLT104</b>	Large	100 gloves/bx	10 bxs/cs
	<b>SCLT105</b>	X-Large	100 gloves/bx	10 bxs/cs

## Certified Claims



SAMPLES: Available by request - ask your sales rep!

\*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

 Licensed for use in Canada

## PEOPLE, PROTECTING PEOPLE

**IMPORTANT NOTE:** The latest product information is at [www.gotgloves.com/products/](http://www.gotgloves.com/products/). Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global