


semper
med

rose™



Care, Comfort
and a Touch of Color



Hospitals



Dental



Veterinary



Ambulatory
Surgical
Centers



Primary
Care



Beauty &
Tattoo

sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





rose™

Sempermed® rose™ exam nitrile gloves are designed to elevate comfort and protection in any setting, from exam rooms to tattoo shops to day spas. These gloves feature a soothing aloe coating that helps prevent skin irritation, keeping hands moisturized and comfortable during extended use. The unique formulation of these gloves enhances protection while maintaining exceptional tactile sensitivity, allowing for precise handling without compromising comfort.

- Aloe & vitamin E inner coating**
 Helps moisturize and care for dry hands, supporting skin comfort during extended glove use.
- Fun unique color**
 Distinctive pink color creates a friendly, professional appearance and helps bring a positive touch to clinical and wellness environments.
- Protects users and patients from the transmission of infections**
 Designed to provide protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

Technical data






TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Pink	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHelf LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.9mils	Palm 2.8mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	542% (avg.)	
BARRIER PERFORMANCE	AQL < 1.5 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN ISO 21420	

Packaging

DIMENSIONS

Dispenser: 245L x 125W x 74H mm
Carton: 385L x 255W x 255H mm

Reorder numbers / Glove Counts*

	SMNR201	X-Small	200 gloves/bx	10 bxs/cs
	SMNR202	Small	200 gloves/bx	10 bxs/cs
	SMNR203	Medium	200 gloves/bx	10 bxs/cs
	SMNR204	Large	200 gloves/bx	10 bxs/cs
	SMNR205	X-Large	180 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global