

semper
med

safe™



**Strong and
Durable**



Hospitals



Lab &
Pharma



Veterinary



First
Responder



Primary
Care



Beauty &
Tattoo

sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





safe™

Sempermed® safe™ exam nitrile gloves offer enhanced safety with increased durability, providing superior chemical resistance, making them a reliable choice in environments with exposure to harsh chemicals. The right combination of protection, comfort, and allergen-free materials, these gloves are an excellent choice for professionals across multiple industries.

- Excellent chemical resistance**
 Against a wide array of chemicals, offers the user enhanced protection and safety.
- Tested against 56 chemotherapy and hazardous drugs, fentanyl, and gastric acid**
- The textured surface on the glove's fingertips offers a strong grip in wet or dry environments.
- Protects users and patients from the transmission of infections**
 Designed to provide a high level of protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.
- ESD tested and approved**
 Safe for applications handling sensitive electronics according to ANSI/ESD STM 15.1

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Blue	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 5.1mils	Palm 3.5mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	564% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, ISO 18889, ASTM D6978	

Packaging

DIMENSIONS
 Dispenser: 245L x 123W x 54H mm
 Carton: 283L x 257W x 261H mm

Reorder numbers / Glove Counts*

	BLNF102	Small	100 gloves/bx	10 bxs/cs
	BLNF103	Medium	100 gloves/bx	10 bxs/cs
	BLNF104	Large	100 gloves/bx	10 bxs/cs
	BLNF105	X-Large	100 gloves/bx	10 bxs/cs
	BLNF106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



avg. 10° D under worst-case test conditions.

SAMPLES: Available by request - ask your sales rep!
 *glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

PEOPLE, PROTECTING PEOPLE