

S semper
med

velvet XP™



Optimal Fit, Extra Protection



Hospitals



Ambulatory
Surgical
Centers



Veterinary



Long Term
Care



Primary
Care



Lab &
Pharma

sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





velvet XP™

Sempermed® velvet XP™ delivers the right combination of fit, feel, tactile sensitivity, strength, and value. This glove is a popular choice for users in a variety of environments and suitable for a wide range of applications. The soft nitrile formulation allows the glove to conform to the hand like a second skin, combining optimal fit and comfort with unrestricted hand motion.

● **Soft as velvet**

The soft formulation of the glove offers enhanced comfort and fit.

● **Your daily choice**

Suitable for a wide variety of medical applications where protection is critical.

● **Tested against 56 chemotherapy and hazardous drugs, fentanyl, and gastric acid**

Meets USP 800 requirements for PPE compliance for chemotherapy.

● **Protects users and patients from the transmission of infections**

Designed to provide a high level barrier of protection against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Violet Blue	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.9mils	Palm 2.8mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	575% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, D6978	

Packaging

DIMENSIONS	250ct	Dispenser: 245L x 128W x 83H mm Carton: 430L x 266W x 257H mm
	200ct	Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm

Reorder numbers / Glove Counts*

	SMTN251 TTNF201	X-Small	250 gloves/bx 200 gloves/bx	10 bxs/cs
	SMTN252 TTNF202	Small	250 gloves/bx 200 gloves/bx	10 bxs/cs
	SMTN253 TTNF203	Medium	250 gloves/bx 200 gloves/bx	10 bxs/cs
	SMTN254 TTNF204	Large	250 gloves/bx 200 gloves/bx	10 bxs/cs
	SMTN255 TTNF205	X-Large	225 gloves/bx 180 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

Velvet 250ct. licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global