

S semper
guard

vinyl comfort™



**All Comfort,
No Latex**



**General
Purpose**



**Janitorial &
Sanitation**



**Safety &
Protection**



**Beauty &
Tattoo**



**Food &
Beverage**

*semperguard is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





vinyl comfort™

Semperguard® vinyl comfort™ gloves are designed to meet the needs of those seeking light duty hand protection at an exceptional value. These polymer coated gloves are manufactured with polyvinyl chloride (PVC) which eliminates exposure to latex and reduces the potential risk of a latex protein reaction. They have a smooth surface and are 100% DEHP and DINP free.

- Compliant with Prop 65 and PFAS-Free**
 Produced without forever chemicals.

- Polymer coated**
 Helping to reduce friction and improving consistent fit in demanding environments.

Technical data






TYPE	Non-sterile vinyl examination and disposable protective glove	
GRADE	General Purpose	
COLOR	Clear	
MATERIAL	Polyvinyl Chloride (PVC)	
GLOVE SURFACE	Smooth	
COUNTRY OF ORIGIN	China	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.1mils	Palm 2.4mils
	based on average production value measured in single layer	
TENSILE STRENGTH <i>(before aging)</i>	19MPa (avg.)	
ELONGATION <i>(before aging)</i>	406% (avg.)	
BARRIER PERFORMANCE	AQL < 4.0 (pinholes)	
SOCIAL COMPLIANCE	EcoVadis	
TESTING & CONFORMITY TO STANDARDS	EN ISO 21420	

Packaging

DIMENSIONS

Dispenser: 200L x 110W x 60H mm
Carton: 310L x 230W x 210H mm

Reorder numbers / Glove Counts*

	GSVF102	Small	100 gloves/bx	10 bxs/cs
	GSVF103	Medium	100 gloves/bx	10 bxs/cs
	GSVF104	Large	100 gloves/bx	10 bxs/cs
	GSVF105	X-Large	100 gloves/bx	10 bxs/cs
	GSVF106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!
**glove counts "by weight"*

Gloves are intended for single use. For maximum protection we recommend changing vinyl gloves between tasks.

PEOPLE, PROTECTING PEOPLE

rev 03/2026 US

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global