



vinyl syn-stretch



The Smarter Stretch



Long Term Care



Veterinary



Lab & Pharma



Primary Care

sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





vinyl syn-stretch

Sempermed® vinyl syn-stretch gloves are designed as the ultimate solution for those who are at risk of becoming sensitive to natural rubber latex. Its cream color and new resin are an upgrade from regular vinyl gloves, giving the look and feel of latex without the latex allergens. It offers exceptional strength and protection, yet maintains maximum tactile sensitivity.

- Compliant with Prop 65 and PFAS-Free**
 Produced without forever chemicals.

- Polymer coated**
 Helping to reduce friction and improving consistent fit in demanding environments.

Technical data

TYPE	Non-sterile vinyl examination and disposable protective glove	
GRADE	Medical	
COLOR	Yellow	
MATERIAL	Polyvinyl Chloride (PVC)	
GLOVE SURFACE	Smooth	
COUNTRY OF ORIGIN	China	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 4.1mils	Palm 3.3mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	23MPa (avg.)	
ELONGATION <i>(before aging)</i>	469% (avg.)	
BARRIER PERFORMANCE	AQL < 1.5 (pinholes)	
SOCIAL COMPLIANCE	EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D5250, EN ISO 21420	

Packaging

DIMENSIONS

Dispenser: 240L x 122W x 60H mm
Carton: 315L x 252W x 256H mm

Reorder numbers / Glove Counts*

	SCVNP102	Small	100 gloves/bx	10 bxs/cs
	SCVNP103	Medium	100 gloves/bx	10 bxs/cs
	SCVNP104	Large	100 gloves/bx	10 bxs/cs
	SCVNP105	X-Large	100 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!
*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing vinyl gloves between tasks.

Licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global