



Semper
med

syntegra CR



70 years of proven performance

Healthcare professional have come to rely on the non-latex comfort and function of polychloroprene.



Easy to don



Optimum protection



Skin-friendly

syntegra CR



Easy to don

Sempermed® syntegra CR is quick and easy to don thanks to an inner coating that considerably reduces friction resistance when putting on a glove. Our sempermed surgical gloves allow for the easy donning of two pairs of gloves. **quick don™**



Optimum protection

The micro-rough surface allows for the confident handling of surgical instruments even in wet environments. The innovative polyisoprene material provides tactile sensitivity comparable to natural rubber latex.



Skin-friendly

Innovative accelerators such as special dithiocarbamates and xanthogenates are used to produce the sempermed® syntegra CR. This reduces the risk of developing type IV allergies.

Technical data

Type	sterile surgical glove for single use, powder free, with synthetic inner coating	
Color	natural	
Size / Overall length as per EN 455-2	6 and 6½ / median 300 mm 7 and 7½ / median 300 mm 8 and 8½ / median 300 mm	
Wall thickness based on avg. production value measured in single layer	Finger .15 mm 6.2 mils	Palm .12 mm 5.0 mils
Barrier performance	AQL 0.65	
Material	polychloroprene	
Glove shape	fully anatomical with curved fingers and rolled rim	
Tensile strength (after aging)	17 MPa min.	
Shelf life in original packaging if stored according to product specification	3 years	
Sterilization	gamma STERILE R	
Labeling & conformity to standards	EN 455-1/-2/-3/-4, EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, ISO 374-5, EN 421, EN 556, ISO 2859-1, ISO 15223-1, ASTM F1671/F1671M	

Purpose

Single-use medical device class IIa
according to MD Directive 93/42/EEC



Single-use protective glove –
PPE category III
(protection against chemical
substances for limited time)*
according to PPE Regulation (EU) 2016/425



Single-use medical device FDA Class I



Packaging

- Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in peel pack.
- In dispenser carton with sterilization indicator: 40 pairs
- In transport carton with sterilization indicator: 240 pairs

Re-order numbers

size 6	SIR600
size 6½	SIR650
size 7	SIR700
size 7½	SIR750
size 8	SIR800
size 8½	SIR850