

velvet XP white



**Comfort and hygiene for
everyday clinical care**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global Brand.**

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet XP white

sempermed® velvet XP white is a powder-free nitrile examination glove designed for reliable barrier protection and precise handling in clinical environments. The glove combines a comfortable fit with dependable tactile sensitivity and a clean white appearance for professional use. An additional variant with extended cuff length (velvet XP white+) provides increased coverage for situations requiring enhanced splash protection.

- White color**
 For a professional and hygienic appearance.
- Textured fingertips**
 Offer secure grip at all times, even in damp conditions.
- Elastic nitrile formulation**
 For comfortable fit and tactile control.
- Extended cuff option (velvet XP white+)**
 For increased coverage and additional splash protection.

Technical data

| | | |
|---|--|--|
| Type | non-sterile, powder-free examination glove for single use | |
| Color | white | |
| Size / Overall length as per EN 455-2 | XS, S, M, L, XL | <i>standard</i> : median ≥ 240 mm <i>velvet XP+ white</i> : median ≥ 265 mm |
| Wall thickness measured in single layer | Finger min. 0.07 mm typical 0.10 mm | Palm min. 0.06 mm typical 0.07 mm |
| Barrier performance | AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement) | |
| Material | Nitrile Butadiene Rubber (NBR) not made with natural rubber latex | |
| Glove surface | finger textured | |
| Glove shape | ambidextrous with rolled rim | |
| Force at break as per EN 455-2 | median ≥ 6 N | |
| Shelf life | 5 years in original packaging if stored according to product specification | |
| Labeling & conformity to standards | EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M | |

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions 150 Dispenser: 240L x 122W x 65H mm
Carton: 340L x 252W x 250H mm

200 Dispenser: 245L x 125W x 74H mm
Carton: 380L x 258W x 255H mm

Available in dispenser box XS, S, M, L 150 / XL 135 pcs
in dispenser box XS, S, M, L 200 / XL 180 pcs
10 dispenser boxes per shipping carton

Article numbers

| | velvet XP white | | velvet XP+ white | |
|-----------|-----------------|------------|------------------|------------|
| XS | 200 pcs | 3000016021 | 150 pcs | 3000016532 |
| S | 200 pcs | 3000016022 | 150 pcs | 3000016533 |
| M | 200 pcs | 3000016023 | 150 pcs | 3000016534 |
| L | 200 pcs | 3000016024 | 150 pcs | 3000016535 |
| XL | 180 pcs | 3000016025 | 135 pcs | 3000016536 |