

**S** semper  
med

# safe+



## Extra long. Extra safe.



Hospital



Primary  
care



Veterinary

Sempermed is a HARPS Global Brand.\*

**CONTACT US!**  
[sempermed@harpsglobal.com](mailto:sempermed@harpsglobal.com)  
[www.sempermed.com](http://www.sempermed.com)

**HARPS Investment Asia Pte. Ltd.**  
1 Harbourfront Avenue, #14-07 Keppel Bay Tower  
Singapore 098632, Tel.: +65 67976560

**HARPS Europe GmbH**  
Wiedner Guertel 9-13,  
1100 Vienna, Austria  
Tel.: +43 1 9346348





# safe+

The sempermed® safe+ provides reliable hand protection when extra barrier and splash resistance are needed. Thanks to the glove's length, wrists and forearms are better protected, which makes it a suitable choice for glove users who work with risk patients or handle substances that can be harmful to the skin.

## Enhanced safety

Additional splash protection thanks to the extra-long cuff.

## Chemical resistance

Barrier protection against certain liquids and chemicals.

## Work with confidence

Finger texturing and good tactile sensitivity.

## Suitable for sensitive skin

Nitrile material helps to prevent type-1 allergies.

## Technical data

<b>Type</b>	non-sterile, powder-free examination glove for single use	
<b>Color</b>	blue	
<b>Size / Overall length</b> as per EN 455-2	XS,S, M, L, XL / median $\geq 295$ mm $\pm 5$ mm	
<b>Wall thickness</b> measured in single layer	Finger min. 0.12 mm typical 0.14 mm	Palm min. 0.08 mm typical 0.09 mm
<b>Barrier performance</b>	AQL $\leq 1.0$ (exceeding the EN 455-1 barrier requirement)	
<b>Material</b>	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
<b>Glove surface</b>	finger textured	
<b>Glove shape</b>	flat with rolled rim	
<b>Force at break</b> as per EN 455-2	typical $\geq 10$ N	
<b>Shelf life</b>	5 years in original packaging if stored according to product specification	
<b>Labeling &amp; conformity to standards</b>	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 421, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

## Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)\* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



## Packaging

<b>Dimensions</b>	Dispenser: 280L x 125W x 60H mm Carton: 310L x 260W x 90H mm
<b>Available</b>	in dispenser box XS, S, M, L 100 / XL 90 pcs 10 dispenser boxes per shipping carton

## Article numbers

<b>XS</b>	100 pcs	3000016026
<b>S</b>	100 pcs	3000016027
<b>M</b>	100 pcs	3000016028
<b>L</b>	100 pcs	3000016029
<b>XL</b>	90 pcs	3000016030