

velvet XP



**Fits like a
second skin**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global brand**

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet XP

sempermed® velvet XP is a powder-free nitrile examination glove designed for everyday use across a wide range of applications. It combines a soft and flexible material with reliable barrier protection and tactile control, supporting comfortable wear that fits like a second skin during routine procedures. Compared to the standard velvet version, it offers slightly increased material strength for users who prefer an added level of protection in daily practice.

■ **Soft and flexible material**

For a comfortable fit and tactile control.

■ **Xtra Protection (XP)**

Slightly increased material strength for users preferring an added level of protection.

■ **Textured fingertips**

For reliable grip, even in damp conditions.

■ **Protects users and patients from the transmission of infections**

Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	violet blue	
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.07 mm typical 0.10 mm	Palm min. 0.06 mm typical 0.07 mm
Barrier performance	AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	flat with rolled rim	
Force at break as per EN 455-2	typical ≥ 7 N	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions	100 pcs	Dispenser: 245L x 125W x 42H mm Carton: 260L x 223W x 255H mm
	200 pcs	Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs in dispenser box XS, S, M, L 200 / XL 180 pcs 10 dispenser boxes per shipping carton	

Article numbers

XS	100 pcs	3000015971	200 pcs	3000015991
S	100 pcs	3000015972	200 pcs	3000015992
M	100 pcs	3000015973	200 pcs	3000015993
L	100 pcs	3000015974	200 pcs	3000015994
XL	90 pcs	3000015975	180 pcs	3000015995

*Sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd., an affiliate of HARPS Global Pte. Ltd.

**Please find below chemical resistance tables for our different glove brands: <https://www.sempermed.com/chemical-resistance-brands>

IMPORTANT NOTES: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **HARPS Global is not liable for incorrect use of the gloves nor will HARPS Global be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classifications set out in this datasheet are up to date as of the date of publication. Where typical technical details are listed, they are lot-specific values from production and may vary in individual cases. Test reports are available upon request. This disclaimer is subject to mistakes, printing errors and amendments. ©Copyright 2026 HARPS Global