

**S** semper  
med

# skySense



Your daily choice of  
nitrile examination gloves

In line with Regulation (EU) 2022/1031 on the International Procurement Instrument – IPI, skySense can be sourced from Malaysia and China under the same article code.



Hospital



Dental



Primary  
care



Long-term  
care



Veterinary

Sempermed is a HARPS Global brand.\*

\* Sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd.

**CONTACT US!**  
sempermed@harpsglobal.com  
www.sempermed.com

**HARPS Investment Asia Pte. Ltd.**  
9 Straits View, #08-10A Marina One West Tower  
Singapore 018937, Singapore  
Tel.: +65 6274 4861 · Fax: +65 6274 6977

**HARPS Europe GmbH**  
Wiedner Guertel 9-13,  
1100 Vienna, Austria  
Tel.: +43 1 9346348

sempermed® skySense offers a lightweight, comfortable and resilient solution for various applications.

The lightweight nitrile examination and protective gloves are available in three colors: blue, violet blue and white. The gloves are not only strong, but also designed to support unrestricted movement and ease of use. They deliver functionality and comfort at great value.

- **Your daily choice:** suitable for a wide variety of medical applications where protection is essential.
- **Tested for protection against chemotherapy drugs.\*\***
- **Protects users and patients from the transmission of infections.** Designed to provide a reliable barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination in clinical settings.

## Technical data

<b>Type</b>	non-sterile, powder-free examination glove for single use
<b>Color</b>	blue / violet blue / white
<b>Size / Overall length</b> as per EN 455-2	XS, S, M, L, XL / median 240 mm
<b>Barrier performance</b> as per EN 455-1	AQL 1.5
<b>Material</b>	Nitrile Butadiene Rubber (NBR)
<b>Glove surface</b>	finger textured
<b>Glove shape</b>	flat with rolled rim, ambidextrous
<b>Force at break</b> as per EN 455-2	median ≥ 6 N
<b>Shelf life</b> in original packaging if stored according to product specification	3 years
<b>Labeling &amp; conformity to standards</b>	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

## Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III (protection against chemical substances for limited time)\*\* as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



## Packaging

available in 100 and 200 pieces per dispenser box  
in dispenser box XS, S, M, L 100 / XL 90 pieces  
in dispenser box XS, S, M, L 200 / XL 180 pieces  
10 dispenser boxes per shipping carton

## Dimensions & Loading Qtys.

<b>DIMENSIONS</b>	Dispenser 245L x 125W x 42H mm
100 pcs.	Carton 260L x 223W x 255H mm
200 pcs.	Dispenser 245L x 125W x 74H mm
	Carton 380L x 258W x 255H mm
<b>LOADING QUANTITIES</b>	20ft 2,033 / 1,187 cartons
	40ft 4,175 / 2,483 cartons
	100/200 pcs. 40ft HC 4,639 / 2,760 cartons

DM1/17 Version 2 01/2026 EN

Article Numbers	skySense	skySense VB	skySense white
<b>XS</b> 100pc 3000016142	200pc 3000015792	100pc 3000016463 200pc 3000015800	200pc 3000016489
<b>S</b> 100pc 3000016143	200pc 3000015793	100pc 3000016464 200pc 3000015821	200pc 3000016490
<b>M</b> 100pc 3000016144	200pc 3000015794	100pc 3000016465 200pc 3000015822	200pc 3000016491
<b>L</b> 100pc 3000016145	200pc 3000015795	100pc 3000016466 200pc 3000015823	200pc 3000016492
<b>XL</b> 90pc 3000016146	180pc 3000015796	90pc 3000016467 180pc 3000015824	180pc 3000016493

\*\*According to the chemical resistance tables for our different glove brands: <https://www.sempermed.com/chemical-resistance-brands>

**IMPORTANT NOTE:** The latest product information is available at [www.sempermed.com](http://www.sempermed.com) under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **HARPS Global is not liable for incorrect use of the gloves nor will HARPS Global be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **HARPS Global bears no liability for allergic reactions.** \*Copyright 2025 HARPS Global

# skySense VB

sempermed® skySense offers a lightweight, comfortable and resilient solution for various applications. The lightweight examination and protective gloves are available in three colors: blue, violet and white. The gloves are not only strong through are designed to support unrestricted and movement and ease of use, they deliver functionality and comfort at a great value.

- **Your daily choice:** suitable for a wide variety of medical applications where protection is critical.
- **Tested for protection against 32 chemotherapy drugs.**
- **Protects users and patients from the transmission of infections.** Designed to provide a high level of protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

## Technical data

<b>Type</b>	non-sterile, powder-free examination glove for single use
<b>Color</b>	violet blue
<b>Size / Overall length</b> as per EN 455-2	XS, S, M, L, XL / median 240 mm
<b>Wall thickness</b> measured in single layer in the palm area	min. 0.05 mm typical 0.06 mm
<b>Barrier performance</b> as per EN 455-1	AQL 1.5
<b>Material</b>	Nitrile Butadiene Rubber (NBR)
<b>Glove surface</b>	finger textured
<b>Glove shape</b>	flat with rolled rim, ambidextrous
<b>Force at break</b> as per EN 455-2	median ≥ 6 N
<b>Shelf life</b> in original packaging if stored according to product specification	3 years
<b>Labeling &amp; conformity to standards</b>	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

## Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)\* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



## Packaging

available in 100 and 200 pieces per dispenser box  
in dispenser box XS, S, M, L 100 / XL 90 pieces  
in dispenser box XS, S, M, L 200 / XL 180 pieces  
10 dispenser boxes per shipping carton

## Dimensions & Loading Qtys.

<b>DIMENSIONS</b>	100 pcs.	Dispenser 245L x 125W x 42H mm Carton 260L x 223W x 255H mm
	200 pcs.	Dispenser 245L x 125W x 74H mm Carton 380L x 258W x 255H mm
<b>LOADING QUANTITIES</b>	20ft	2,033 / 1,187 cartons
	40ft	4,175 / 2,483 cartons
	100/200 pcs.	40ft HC 4,639 / 2,760 cartons

## Article numbers

<b>XS</b>	100pc	3000016463	200pc	3000015800
<b>S</b>	100pc	3000016464	200pc	3000015821
<b>M</b>	100pc	3000016465	200pc	3000015822
<b>L</b>	100pc	3000016466	200pc	3000015823
<b>XL</b>	90pc	3000016467	180pc	3000015824

DM1/17 Version 1 07/2025 EN

# skySense white

sempermed® skySense is a lightweight examination and protective glove. These gloves are flexible and resilient, providing a lightweight yet strong solution for various applications. Designed for unrestricted hand movement and ease of use, combining functionality and comfort all at a great value.

- Your daily choice:** suitable for a wide variety of medical applications where protection is critical.
- Tested for protection against 32 chemotherapy drugs.**
- Protects users and patients from the transmission of infections.** Designed to provide a high level of protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

## Technical data

<b>Type</b>	non-sterile, powder-free examination glove for single use
<b>Color</b>	violet blue
<b>Size / Overall length</b> as per EN 455-2	XS, S, M, L, XL / median 240 mm
<b>Wall thickness</b> measured in single layer in the palm area	min. 0.05 mm typical 0.06 mm
<b>Barrier performance</b> as per EN 455-1	AQL 1.5
<b>Material</b>	Nitrile Butadiene Rubber (NBR)
<b>Glove surface</b>	finger textured
<b>Glove shape</b>	flat with rolled rim, ambidextrous
<b>Force at break</b> as per EN 455-2	median ≥ 6 N
<b>Shelf life</b> in original packaging if stored according to product specification	3 years
<b>Labeling &amp; conformity to standards</b>	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

## Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)\* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



## Packaging

available in 200 pieces per dispenser box  
in dispenser box XS, S, M, L 200 / XL 180 pieces  
10 dispenser boxes per shipping carton

## Dimensions & Loading Qtys.

<b>DIMENSIONS</b>	Dispenser 245L x 125W x 74H mm Carton 380L x 258W x 255H mm
<b>LOADING QUANTITIES</b>	20ft 1,187 cartons 40ft 2,483 cartons 40ft HC 2,760 cartons

## Article numbers

<b>XS</b>	200pc	3000016489
<b>S</b>	200pc	3000016490
<b>M</b>	200pc	3000016491
<b>L</b>	200pc	3000016492
<b>XL</b>	180pc	3000016493

DM1/17 Version 1 07/2025 EN